

Access And Comfort  
Are Easily Within Reach.  
The Performer Chair By A-dec.



**a dec**<sup>®</sup>

# More Choices, Lower Price, Greater Value.



*Performer's thin-profile contoured backrest upholstery offers more leg room and is wide enough to comfortably support a variety of shoulder widths.*



*Convenient for hygiene and pediatric dentistry, Performer's optional narrow backrest improves access to the patient—especially while working direct.*



*To accommodate a wider variety of shoulder widths, yet still allow for excellent direct access, choose the narrow backrest with slings option.*



*Performer's chair base features A-dec's proven hydraulic cantilever lift system for smooth and quiet operation.*



*A programmable chair touch pad is available when you order the Performer chair with an integrated Performer delivery system.*

The Performer® III chair offers you more flexible choices, better access, and improved comfort—all at a price you can afford. And since Performer is designed, built, and tested by A-dec, you can also expect to get the most out of every dollar you invest.

## FLEXIBLE CHOICES FOR COMFORT AND ACCESS

Performer's flexible options are designed to save you money. They're also designed to help you do better dentistry. The chair's thin-profile upholstery choices, for instance, allow you to lower the patient more, so you can bring your arms down to a more comfortable working position. For treating smaller patients and working direct, you can choose Performer's narrow backrest, which helps enhance your comfort by allowing you to move closer to the patient's shoulders during procedures.

Since the patient's comfort is just as important as yours, Performer's backrest, seat, and toeboard are designed to help cradle the patient and reduce pressure points. Performer is also available with slings, so you can choose the narrow backrest and still comfortably support most patients.

## BACKED BY A-DEC INNOVATION AND QUALITY

With three backrest choices, two different mounting points for the delivery system, and a number of attachments, Performer offers you more innovative choices than any other chair in its price range. But most important, Performer is built to last. The chair includes the same hydraulic cantilever systems, for example, as Cascade® and Decade®—two of the most popular chairs in dentistry. This proven system features lift cylinders built entirely by A-dec, so Performer operates quietly, smoothly, and dependably. Consider the rest of Performer's quality A-dec components, and you'll begin to discover a chair that's not only lower in cost, but also much higher in value.

For a closer look at the Performer chair, contact your local authorized A-dec dealer, or call **1-800-547-1883** for a list of dealers in your area.

Product Features	Standard	Optional
Thin-Profile Backrest	1	
Narrow Backrest Upholstery		m
Slings		m
Double-Articulating Headrest	1	
Pivoting Armrests		m
Programmable Button Footswitch	1	
Programmable Chair Touch Pad (on delivery systems)		m
Cast-Aluminum Baseplate	1	
Hydraulic-Lift Cantilever System	1	
Accessory Power Cord		m
Multiple Mounting Points (for system, cuspidor, light, and other options)	1	



A-dec Inc.  
2601 Crestview Drive  
Newberg, OR 97132 USA  
Phone: 1-800-547-1883, Fax: 1-503-538-0276

Outside USA/Canada, contact:  
A-dec International Inc.  
Phone: 1-503-538-7478, Fax: 1-503-538-5911  
©A-dec Inc., 2001. All rights reserved. 85.6067.00/Dyna/15M/11-01

*The colors and products shown in this sales sheet may vary from the actual materials and are subject to change. Please see your authorized A-dec dealer for accurate color samples and the most current product information.*