



Dentaquip

WATER DISTILLERS



- Flash vaporization
- Condensation
- 99% Pure water

SciCan.Aquastat

A new wave breakthrough provides 99% pure, fresh water. Protect your patients and your Statim Cassette Autoclave® by using only pure, fresh water in your practice. SciCan's AQUASTAT distiller allows you to forget about buying expensive bottled water or risking damage to your equipment from less than pure distilled water. A patented distillation technology kills harmful bacteria and viruses by flash vaporizing water and then reconstituting it through condensation. The result is water that's 99% pure. So pure, you could even drink it!

AQUASTAT purifies soft or hard water from any source including those that are spring fed or municipally piped. Simply plug it into any outlet, fill with tap water and switch on. Safely and effectively it will remove lead, copper, aluminum, arsenic, iron, nitrates, herbicides, pesticides, PCBs and more! With its energy efficient technology, AQUASTAT purifies 3.78 litres (1 gallon) of water in just 4.25 hours. Combine this fast process with a compact, light-weight design that doesn't require any special plumbing or wiring - and you'll have all the fresh water.you.need.

SPECIFICATIONS

| | |
|---------------------|--|
| Weight: | 4.5 kg (10 lb.) |
| Dimensions : | 40.6 x 24.5 x 38.4 cm (L x W x H) (16 x 9 5/8 x 15 1/8") |

.....

Dentaquip Limited
 1 Altona Road // Lisburn // N. Ireland// BT27 5QB
 t: 028 9260 1000 // f: 028 9267 6550
 Reg. No. NI 9638

