




**KaVo SUPERtorque.**  
The perfect angle.



KaVo. Dental Excellence.




## KaVo SUPERtorque.

### The perfect angle. The perfect turbine:

The angle is the key. With the new SUPERtorque Turbine 660 I can access even hard-to-reach areas. I have plenty of room to manoeuvre – for whatever the situation or treatment requires. SUPERtorque 660 also comes with a wealth of amazing features such as a firm grip, the durability you expect from a KaVo product and handling comfort that is second to none. The quality you associate with KaVo.

Simply perfect. The new SUPERtorque 660 turbines are the next generation of turbines, combining world-renowned KaVo quality with the latest in KaVo innovation. With its KaVo-patented angle combination (100° head angle and 18° knee angle), the new SUPERtorque turbine range is ideal for precision work particularly in the less accessible molar region of the mouth.



# KaVo SUPERtorque.


## Total silence. Full power.

Torque is a measure of 'force exerted in rotation'. And the name certainly fits here. The SUPERtorque 660 is incredibly well balanced and it is this balancing act between minimalism and high performance, between quiet operation and extreme power that goes somewhere towards explaining the allure of KaVo turbines. I can feel it. The SUPERtorque 660 brings a sense of spontaneous enjoyment to my day-to-day work.

The new SUPERtorque 660 turbines operate at 18 Watt, offering you top performance and maximum power in one. And their vibration-free operation enables safe, precision work in any situation. With 25,000 LUX, the glass rod fibre-optic bundle perfectly illuminates even the least accessible preparation areas. SUPERtorque 660 turbines also have a wealth of excellent standard features that have long been exclusive to KaVo.



**MULTiflex coupling**  
Complete compatibility. SUPERtorque 660 is compatible with all couplings of the KaVo MULTiflex system, which means that it can be used with any treatment unit.



**Hygiene**  
As clean as a whistle. KaVo instruments can be thermally disinfected and sterilised in an autoclave at up to 135 °C.



**Service life**  
On track for the entire practice. KaVo original ceramic ball bearings guarantee a much longer than average service life and excellent, vibration-free operation.





**Cooling**

The temperature is always right. The four-port spray emits a fine mist of spray, ensuring that the instrument is always optimally cooled – right to the bur tip.

**Light**

A perfect view. With 25,000 LUX, the glass rod fibre-optic bundle of KaVo SUPERtorque 660 B illuminates the preparation area perfectly – with not a shadow to be seen.

**Ease of use**

Reliable and handy; with a retention power of 32 N, the KaVo pushbutton chuck system holds the bur securely in place, facilitating complete ease of use.

# KaVo SUPERtorque 660 B / C.

## Just one decision. A galaxy of benefits.

I already have experience of working with KaVo turbines. I want to use KaVo turbines in my work. I will always use KaVo turbines in my work. Whatever experience you have of KaVo turbines and whatever you ask of them, one thing is certain: Choosing SUPERtorque 660 means choosing KaVo over the long term. With all the benefits - or to put it another way - with all the advantages.

### **SUPERtorque 660 B**

*REF 1.002.2200*

### **SUPERtorque 660 C**

*REF 1.002.2100*

- Glass rod fibre-optic bundle (25,000 LUX)\*
- Ideal angle combination (100° head/18° knee)
- Quiet, vibration-free operation (62 dB(A))
- 18 Watt power
- KaVo original ceramic ball bearings
- KaVo pushbutton chuck system (32 N retention power)
- 4 port spray for optimal cooling
- KaVo MULTIflex coupling system
- Can be thermally disinfected
- Can be sterilised in an autoclave at up to 135 °C

\*only available with SUPERtorque 660 B



## KaVo QUATTROcare, KaVo SPRAYrotor, KaVo Spray.

### Enjoy top quality – and care that keeps it that way.

Every KaVo instrument is of the highest quality. And I want it to stay that way. To prevent KaVo instruments from losing any of their lustre or value, I care for them using a special, KaVo-designed care programme. This programme includes KaVo QUATTROcare, KaVo SPRAYrotor and KAVO spray. Durability does not depend on quality alone; it also depends on my approach to conserving value.

To accompany its extremely high quality instruments, KaVo has a range of equally high quality care products. These products have been specifically designed to suit the care requirements of each individual instrument; they are of the highest quality and guarantee a long service life for your instruments.

**KaVo QUATTROcare 2104 A**  
Value conservation par excellence. The innovative care system automatically cleans and maintains up to four instruments at a time.  
REF 0.411.7800

**KaVo SPRAYrotor 3 N**  
Intensive cleaning. Rotates to reach even the deepest, narrowest recesses.  
REF 0.411.7710

**KaVo Spray**  
Efficient cleaning. 500ml for manual maintenance of all KaVo straight and contra-angle handpieces, heads and turbines.  
REF 0.411.9640





D-88400 Biberach/Riß  
Telefon: +49 73 51 56-0  
Fax: +49 73 51 56-14 88  
Internet: [www.kavo.com](http://www.kavo.com)