


Note: Please note all promotion prices shown have already had the trade in discounts deducted. The prices will not be further reduced.

SUPERtorque 660B with fibre optics and 660C without fibre optics

Efficient Proven Turbine Technology

The best known and best loved turbine available with a long running legacy for quality and reliability

- Push button chuck
- 100 degree working angle
- 4 port spray
- 25.000 Lux fibre optics for excellent illumination - 660B only




Trade in your old turbine and get £125 off your KaVo 660B

The KaVo SUPERtorque 660B for only £660 (RRP £945)
The KaVo SUPERtorque 660C for only £565 (RRP £832)
Pound for pound no other turbine is better value for money!

GENTLesilence 8000 B with fibre optics and 8000 C without fibre optics

- Incredibly quiet at only 57 dB(A) with dampened frequency range
- Increased power at 19 watts
- Small dimensions (head 13.4mm) for improved access and added patient comfort
- Plasmatec surface for improved comfort and hygiene
- Absolutely illuminating with 25.000 Lux fibre optics – 8000B only
- Ideal head angle giving more space in confined preparation areas
- Ceramic ball bearings
- Washable in Thermo disinfectant



Trade in your old turbine and get £150 off your KaVo 8000B

Promotions: The Ultimate Turbines at massively reduced prices:

8000B for only £815 (RRP £1,205)
8000C for only £725 (RRP £1,060)

Quattrocare 2104

Automated maintenance system for cleaning and care of handpieces

Designed for trouble free cleaning, purging and lubrication of up to four handpieces simultaneously within 1 minute

A free SUPERtorque 660B with every QUATTROcare, worth £945

Promotion Price: £1995 – A small price to pay to look after your handpiece investment - savings of over £1,000



SONICflex 2003 Air Scaler without light

- Frequency: 6000 Hz
- 3 performances levels 120µm, 160µm, 240µm
- Direction of movement: elliptical
- Very quiet
- Includes 3 tips

Promotion: **Buy a SONICflex 2003 for only £882 (RRP £1,290)**



**Trade in your old
Air Scaler and get
£100 off your
KaVo SONICflex**

All prices exclude VAT

All prices are shown as the final promotion price after deduction of the trade in value.